

SUGGESTIONS TO HELP WITH THE PATCH IRRITATION AND REMOVAL

1. Prior to applying the patch use either Vitamin E - pierce the capsule and use the gel - or a liquid anti-acid (Mylanta or Milk of Magnesia) and apply a thin coat around the eye. Let dry and then put on the patch.
2. Prior to putting patch, put the patch on your arm or pant leg so some of the adhesive gets "contaminated" and is less sticky.
3. If there are certain areas that have the sores, cut the adhesive around these areas if possible so they don't stick there but do everywhere else.
4. Alternate turning the patch small side in toward nose, then next day large size toward nose if this keeps the adhesive from the area with the sores.
5. When removing the patch, we always coat the outside of the patch with Vaseline or baby oil, Aquaphor or diaper ointment about 15 minutes before we remove it to loosen the adhesive abit.

These ideas have been provided to us by parents and doctors

Fresnel Prism and Lens does have available for purchase a product that provides an adhesive barrier that last for up to 3 days before needing to be reapplied. Cavilon allows the skin to heal while remaining patched.

*MYI™ Occlusion Therapy Patches
By the Fresnel Prism & Lens Co.
10810 Nesbitt Ave S.
Eden Prairie, MN 55437 USA*

*p: 952-496-0432 or 1-800-544-4760 US & Canada
e: info@fresnelprism.com
w: www.myipatches.com*